

TITRATION OF AAV8 RSM VECTOR GENOMES USING QUANTITATIVE REAL TIME PCR

Place 5 μ L of rAAV preparation and controls (PBS and rAAV control) in the 1.5ml eppendorf tube.

Prepare a DNase digestion mix:

	<u>[Final]</u>
5 μ l 10x DNase I buffer	1X
0.5 μ l 100X PF68 (5% w/v)	0.05%
2 μ l DNase I (10u/mL)	400U/mL
37.5 μ l NF H ₂ O	

Add 45 μ L of DNase digestion mix per tube. Incubate 60 minutes at 37°C.

Prepare Samples

Prepare Sample dilution buffer:

	<u>[Final]</u>
1mL 1X GeneAmp PCR Buffer	1X
0.1mL 100X SSS DNA 200 μ g/mL	2 μ g/mL
0.1mL 100X pluronic F68 (5%)	0.05%
8.8mL NF H ₂ O	

In duplicate, prepare a serial dilution of extracted sample and controls: Dilutions 1 and 2: 90 μ L sample dilution buffer + 10 μ L sample (10-fold dilution) Dilutions 3 to 7: 40 μ L sample dilution buffer + 10 μ L sample (5-fold dilution)

Prepare standard plasmid

Prepare a 10-fold serial dilution of RSS rAAV plasmid (pTR-UF11) in sample dilution buffer (10 to 10⁸ copies/20 μ L).

See appendix A

Prepare PCR reactions

Prepare enough mix for 2 reactions of each DNA sample, a standard curve, two control samples (PBS + rAAV control), plus 2 reactions for no template control (H₂O).

<u>Reagent</u>	<u>Vol. per 50ul PCR reaction</u>	<u>[Final]</u>
TaqMan Universal PCR Mix (2X)	25 μ L	1X
Forward Primer (10 μ M)	1 μ L	0.2 μ M
Reverse Primer (10 μ M)	1 μ L	0.2 μ M
Fluorescent Probe (10 μ M)	0.5 μ L	0.1 μ M
Nuclease-free water	2.5 μ L	
Sample DNA	20 μ L	

NB: Lower reaction volumes can be used to accommodate your real-time PCR machine – reagents should be scaled down appropriately.

Place the optical caps or optical plate sealer sheet over the plate. Centrifuge briefly before placing in the machine.

TITRATION OF AAV8 RSM VECTOR GENOMES USING QUANTITATIVE REAL TIME PCR

Cycles:	Hold	95	10:00
	Cycle 40	95	00:15
		60	01:00

Analyzing Data

Refer to the appropriate instrument user guide for instructions on how to analyse your data.

The general process for analyzing the data involves the following procedures:


1. View the amplification plots

2. Set the baseline

During early PCR cycles, the background signal in all wells is used to determine the baseline fluorescence. This is usually cycles 3-12.

3. Set the threshold

The Threshold should be placed in the region of exponential phase. To standardize threshold placement we recommend you adjust the threshold until the 10^8 copies standard reads ~ 13.08 Cts. The replicates should be tight, within 0.5Cts (Threshold Cycle).



Plasmid range analysis


Plot the Cts of your standard plasmid (Y axis) versus the log of your initial quantity (X axis) to generate a standard curve. The standard curve should be linear over the entire range of where you expect your unknowns to fall. Standard curve specifications should be:

TITRATION OF AAV8 RSM VECTOR GENOMES USING QUANTITATIVE REAL TIME PCR

- Correlation Coefficient ≥ 0.99
- Efficiency of amplification 95-105%

Example:

Copy number	ln (copy number)	Ct1	Ct2	Ct average
9,00E+07	18,32	12,8	12,9	12,9
9,00E+06	16,01	15,7	14,7	15,2
9,00E+05	13,71	18,2	18,5	18,4
9,00E+04	11,41	21,7	21,7	21,7
9,00E+03	9,10	25,2	25,3	25,2
9,00E+02	6,80	28,6	28,6	28,6
9,00E+01	4,50	32,0	31,9	31,9
9,00E+00	2,20	35,6	35,7	35,6



Regression line: $y = -a \ln(x) + b$

Negatives controls:

Cts of controls should be above 35 (PBS and H₂O).

Determination of titer:

1-Determine number of copies per PCR reaction for each dilution of test article using the standard curve . This feature is automated in most real time PCR software packages. For those performing the calculations manually, determine the number of copies of each sample by applying the equation of the regression line (plasmid range).

$$\text{Copy number} = e^{[(Ct-b)/a]}$$

2- According to the dilution of the sample and volume of sample loaded, calculate the number of copies per ml (vg/ml). NB: since real time PCR only targets one of the two strands packaged within the AAV capsid, a multiplication factor of 2 is also required.

3- Carry out the average of the titers obtained for each dilution to arrive at the final

4-rAAV control: Check that the titer obtained is included in the interval of data previously obtained.

Appendices

A) Preparation of Standard Plasmids for Real-Time PCR

Materials

- 20% Lysol
- 70% Ethanol
- Plasmid
- Restriction Enzyme with Buffer
- Sample dilution buffer (from the protocol above)
- 15mL polypropylene tube
- Vortex mixer
- Biosafety cabinet with UV capabilities
- QIAquick PCR Purification Kit (Qiagen, Cat# 28104)
- QIAmp DNA Mini Kit (Qiagen, Cat# 51304)
- Spectrophotometer

Procedure

Note: Work inside the biosafety cabinet to avoid potential contamination from the laboratory environment.


- Clean the biosafety cabinet work surface and pipette shafts with 20% Lysol, followed by 70% Ethanol. Irradiate with UV light for 10-15 minutes.
- Linearize the standard plasmid with a restriction enzyme that cuts outside your intended PCR target.
- Purify the linearized plasmid using the QIAquick PCR purification kit or similar if it is less than 10kb in size. For larger plasmids, use the QIAmp DNA mini kit or similar. In both cases, pre-warm the elution buffer to 70°C before use. Also, pre-warm the columns loaded with elution buffer at 70°C for 5 minutes. Carry out two elutions in a minimum volume to maximize recovery and concentration of the plasmid.

Note: Do not carry out a phenol:chloroform: isoamyl purification of the digest as any contamination with phenol will affect the accuracy of the spectrophotometric reading.

- Determine the concentration of the linearized plasmid by spectrophotometry and convert the concentration readout to grams per Liter (g/L).
Note: Make sure the spec reading falls within the linear range of the instrument. Dilute the sample less if necessary. Blank the reading against an equal dilution of elution buffer and use the background correction feature of the instrument (read at 320nm).
- Calculate the formula weight (F.W.) of your standard plasmid: F.W. = Plasmid size (in base pairs) x 662 grams/mole.bp **(For TRUF11, this is 4766400)**
- Calculate the molar concentration (M) of your linearized plasmid. $M = \text{mole/L} = (\text{mass (in grams)}/\text{F.W.})/1\text{L}$.
- Determine copy number per microliter of your linearized plasmid based on molar concentration. 1M is equivalent to about 6.0221415×10^{23} copy numbers.
- Make the first dilution, typically with a final concentration of 1×10^{10} copies per 20 μL (5×10^8 copies/ μL). Use the $C_1V_1=C_2V_2$ relationship, where V_1 is the unknown volume, C_1 is the stock linearized plasmid concentration (copies/ μL), V_2 is the final volume (typically 100 μL) and C_2 is the final concentration (5×10^8 copies/ μL). Add V_1 volume of stock linearized plasmid to enough Sample dilution buffer for 100 μL of solution.

TITRATION OF AAV8 RSM VECTOR GENOMES USING QUANTITATIVE REAL TIME PCR

- Carry out a serial dilution of 1×10^{10} copies/ $20 \mu\text{L}$ of the linearized standard using sample dilution buffer. A volume of $90 \mu\text{L}$ of diluent is used per dilution



B) AAV8 Reference Standard Material Test Record Packet.